GlucoSure[®]Aut[®]Code

Blood Glucose Monitoring System











Autocoding



Alternate site testing



Foil pack available

Result in 6 seconds





PC link

Ë

Ë

GlucoSure AutoCode is an unique autocoding blood glucose monitoring system that automatically selects the best calibration each time a test strip is inserted. Calibration error is often the leading cause of inaccurate readings. GlucoSure Autocode eliminates calibration errors and makes testing simpler. No more code card or button-press code selection is required.

GlucoSure AutoCode offers fast and accurate result in 6 seconds using tiny amount of blood. The meter is small and compact and can fit into any life style. The large LCD screen is designed for easy reading and up to 300 test results with date and time can be stored and recalled for later use. The stored test results can be averaged in 7, 14, and 30 days to better understand blood glucose management. The stored results can also be downloaded into a computer for further analysis.



- Measurement Range 20~600 mg/dL
- Coding Method Automatic
- Quick Testing Time 6 seconds
- Memory 300 test results
- Data Management 7, 14 and 30 days average
- Sample Volume 0.6 µL
- Additional Features PC link, Bluetooth (optional)



ISO 15197 : 2015		
Sample < 100 mg/dL		
±5 mg/dL	±10 mg/dL	±15 mg/dL
125 / 180	174 / 180	180 / 180
69%	97%	100.0%
Sample ≥100 mg/dL		
±5%	±10%	±15%
271/420	397 / 420	419 /420
65%	95%	99.8%



No.7, Li-Hsin Rd. V, Hsinchu Science Park, Hsinchu, Taiwan, R.O.C. TEL: 886-3-5641952 / FAX: 886-3-5678302 / E-Mail: info@apexbio.com / Website: www.apexbio.com